

AeroDawg 

Canine Aerosol Chamber




SMALL

MEDIUM

LARGE

Mask Sizer

Each device is shipped with 2 different sized masks.



AeroDawg* Chamber SMALL,
includes small and medium masks

AeroDawg* Chamber LARGE,
includes medium and large masks

How to Use

The mask selected should cover the nose and mouth. However, if not possible apply mask to the dog's nostrils while gently holding the muzzle closed.



Trudell Medical International

www.aerodawg.com



201092-001 Rev.A. * trade-marks and registered trade-marks of Trudell Medical International.
© TMI 2015. All rights reserved.